

Архангельск (8182)63-90-72	Калининград (4012)72-03-81	Нижний Новгород (831)429-08-12	Смоленск (4812)29-41-54
Астана +7(7172)727-132	Калуга (4842)92-23-67	Новокузнецк (3843)20-46-81	Сочи (862)225-72-31
Белгород (4722)40-23-64	Кемерово (3842)65-04-62	Новосибирск (383)227-86-73	Ставрополь (8652)20-65-13
Брянск (4832)59-03-52	Киров (8332)68-02-04	Орел (4862)44-53-42	Тверь (4822)63-31-35
Владивосток (423)249-28-31	Краснодар (861)203-40-90	Оренбург (3532)37-68-04	Томск (3822)98-41-53
Волгоград (844)278-03-48	Красноярск (391)204-63-61	Пенза (8412)22-31-16	Тула (4872)74-02-29
Вологда (8172)26-41-59	Курск (4712)77-13-04	Пермь (342)205-81-47	Тюмень (3452)66-21-18
Воронеж (473)204-51-73	Липецк (4742)52-20-81	Ростов-на-Дону (863)308-18-15	Ульяновск (8422)24-23-59
Екатеринбург (343)384-55-89	Магнитогорск (3519)55-03-13	Рязань (4912)46-61-64	Уфа (347)229-48-12
Иваново (4932)77-34-06	Москва (495)268-04-70	Самара (846)206-03-16	Челябинск (351)202-03-61
Ижевск (3412)26-03-58	Мурманск (8152)59-64-93	Санкт-Петербург (812)309-46-40	Череповец (8202)49-02-64
Казань (843)206-01-48	Набережные Челны (8552)20-53-41	Саратов (845)249-38-78	Ярославль (4852)69-52-93

Единый адрес для всех регионов: <http://unispectr.nt-rt.ru/> || urc@nt-rt.ru

Миксер MTR22

Технические характеристики

Laboratory Mixers

Multi-Purpose Mixer

The MTR22 Multi-Mix has the capability to perform three modes of mixing, all of which are programmable and can all be used in the same "run". Or, if you just want to rotationally mix, set the unit to rotate. It is very simple to operate and program, runs can be programmed in seconds. You can mix many different sizes of tubes ranging from 1.5 ml to 50 ml with the two interchangeable mixing platforms that come standard with the unit. Tubes are held in-place with unique rubber flexible grippers; that surround and firmly hold tubes of varying diameters.

The MTR22 can mix in the following motions:

1. Rotate 360°, forward or reverse, and forward, then reverse
2. Reciprocating Rocker, back and forth like a traditional rocker
3. Vortex, (cavitate)

Tube Capacity, Two Platforms:

Platform 1: 22 Tubes: 1.5 ml – 15 ml (Installed when new)

Platform 2: Six, 50 ml Tubes and Two 1.5 – 15 ml Tubes
(Included as an accessory)

Easy to use controls to set time, shaking speed, rotation angle, and pauses. Use the Multi-Mixer to mix homogenize samples of all kinds, mix blood samples, urines, injectibles, tissue culture, plasma, hybridization and more. Mix anything that needs to be mixed in a tube. The easy to program unit will be a valued asset in your lab! You can load tubes, program and walk away! Let the mixer do the work for you.



Specifications:

- Capacity:
 - Platform 1: 22 Tubes: 1.5 ml – 15 ml (Installed when new)
 - Platform 2: Six, 50 ml Tubes and Two 1.5 – 15 ml Tubes (Included as an accessory)
- RPM: Variable from 5 – 30 RPM
- Timer: Up to 999 Seconds, max mixing time: 12 hours
- Reciprocating Timer: Up to 999 Seconds
- Rotation Angle in Reciprocation Mode: 15° - 90° (in 15° increments)
- Turning angle vortex mode: 1 to 5
- Vortexing times: 1 to 5
- Maximum timer: 12 hours
- Power: 110 / 220 V
- Dimensions: 13.5" x 7" x 7.5"
- Weight: 6.7 lbs. / 3 Kgs

Item No.	Description
MTR22	Programmable Multi-Mix Rotator, 110/220V



MTR22 shown with 6-place 50ml platform installed.



MTR22 includes two interchangeable platforms.

Syringe and Instrument Tray



Great for storing syringes of various sizes and other accessories.

Specifications:

- Six compartments, each measures 2.25"W x 9"L x 1"D
- Overall dimensions: 16"L x 9.5"W x 1.5"D
- Sold in a package of two trays

Item No.	Description
48400	Syringe & Instrument Tray, Pk/2

Архангельск (8182)63-90-72	Калининград (4012)72-03-81	Нижний Новгород (831)429-08-12	Смоленск (4812)29-41-54
Астана +7(7172)727-132	Калуга (4842)92-23-67	Новокузнецк (3843)20-46-81	Сочи (862)225-72-31
Белгород (4722)40-23-64	Кемерово (3842)65-04-62	Новосибирск (383)227-86-73	Ставрополь (8652)20-65-13
Брянск (4832)59-03-52	Киров (8332)68-02-04	Орел (4862)44-53-42	Тверь (4822)63-31-35
Владивосток (423)249-28-31	Краснодар (861)203-40-90	Оренбург (3532)37-68-04	Томск (3822)98-41-53
Волгоград (844)278-03-48	Красноярск (391)204-63-61	Пенза (8412)22-31-16	Тула (4872)74-02-29
Вологда (8172)26-41-59	Курск (4712)77-13-04	Пермь (342)205-81-47	Тюмень (3452)66-21-18
Воронеж (473)204-51-73	Липецк (4742)52-20-81	Ростов-на-Дону (863)308-18-15	Ульяновск (8422)24-23-59
Екатеринбург (343)384-55-89	Магнитогорск (3519)55-03-13	Рязань (4912)46-61-64	Уфа (347)229-48-12
Иваново (4932)77-34-06	Москва (495)268-04-70	Самара (846)206-03-16	Челябинск (351)202-03-61
Ижевск (3412)26-03-58	Мурманск (8152)59-64-93	Санкт-Петербург (812)309-46-40	Череповец (8202)49-02-64
Казань (843)206-01-48	Набережные Челны (8552)20-53-41	Саратов (845)249-38-78	Ярославль (4852)69-52-93

Единый адрес для всех регионов: <http://unispectr.nt-rt.ru/> || urc@nt-rt.ru